

# FREEDOM 4

## 4-Channel 12-Band Digital Amplifier 100% Digital Processing

### Features and available options:

- Adaptive Feedback Rejection for Improved Stability, Higher Gains, and Larger Venting
- New Layered Noise Reduction Reduces Noise When Speech is not present
- 4-Channel WDRC Compression
- Excellent Performance in Noisy Environments
- Selectable Time Constant Modes
- Look-ahead Detection for Loud Impulse Noises
- Independent Expansion in Each Channel
- Optional Auto T-coil or Manual T-coil
- A variety of tube and dome sizes available
- Available in 5 different case colors  
Beige, Gray, Charcoal, Cocoa Brown, and Black
- 12-Band Equalizer for Fine Tuning at all Frequencies
- Up to four Programmable Memories with Memory Identification Tones
- Selection of Directional Patterns
- Wideband Output Compression limiting
- MPO Limiting Threshold Adjustments for each Channel
- Low Battery Warning Tone Alerts User of Low Battery Level


|                              |          |
|------------------------------|----------|
| Battery Size                 | 10A      |
| SSPL90 Max. (db SPL)         | 120      |
| HF Average OSPL90 (db SPL)   | 110      |
| Peak Gain (db)               | 35       |
| HF Average Full on Gain (db) | 30       |
| Frequency Range (Hz)         | 200-8000 |
| Batter Life (Hours)          | 90       |

### Electronic Digital Processing Option:

- Programmable Directional Patterns --
- Omni Directional                      Cardioid
- Super-cardioid                          Hyper-cardioid



**FREEDOM 4  
OPEN FIT MINI-OTE  
(Over the Ear)**



## Highest Quality at Lowest Prices

Please Contact a Prairie Labs Sales Representative for more details (800) 322-8238